

Human VRK3 protein (Recombinant) (His-Tag) (STJP017695)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 40% Glycerol, 0.15M NaCl, 1mM DTT
Molecular Weight	55.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51231
Gene Symbol	VRK3
UniProt ID	VRK3_HUMAN
Immunogen Region	1-474aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSMISFCPD CGKSIQAQAFK FCPYCGNSLP VEEHVGSQTF VNPHVSSFQG SKRGLNSSFE TSPKKVKWSS TVTSPRLSLF SDGDSSSEED TLSSSSERSKG SGRPPTPKS SPQKTRKSPQ VTRGSPQKTS CSPQKTRQSP QTLKRSRVTT SLEALPTGTV LTDKSGRQWK LKSFQTRDNQ GILYEAAPTS TLTCDSGPQK QKFSLKLD

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human VRK3 protein under reducing conditions Coomassie blue stain