

Human NMRAL1 protein (Recombinant) (His-Tag) (STJP017696)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol, 0.15M NaCl
Molecular Weight	35.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57407
Gene Symbol	NMRAL1
UniProt ID	NMRL1_HUMAN
Immunogen Region	1-299aa
Immunogen Sequence	MVDKKL VVVFGGTGAQ GGSVARTLLE DGTFKVRVVT RNPRKKAKE LRLQGAENVQ GDQDDQVIME LALNGAYATF IVTNYWESCS QEVEVKQGKL LADLARRLGL HYVVYSGLEN IKKLTAGRLA AAHFDGKGEV EEFYFRDIGVP MTSVRLPCYF ENLLSHFLPQ KAPDGKSYLL SLPTGDPMD GMSVSDLGPV VLSLLK

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human NMRAL1 protein 3µg reducing conditions Coomassie blue stain