

## Human RhoG protein (Recombinant) (His-Tag) (STJP017705)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

### PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) 20% Glycerol, 50mM Imidazole, 0.5M NaCl
Molecular Weight	25.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">391</a>
Gene Symbol	<a href="#">RHOG</a>
UniProt ID	<a href="#">RHOG_HUMAN</a>
Immunogen Region	1-188aa
Immunogen Sequence	MQS IKCVVVG DGA VGKTCLLICY TTNAFPKEYI PTVFDNYSAQ SAVDGRVTNL NLWDTAGQEE YDRLRTLSTYP QTNVFCVCF S IASPPSYENV RHKWHPEVCH HCPDVPILLV GTKKDLRAQP DTLRRLKEQG QAPITPQQGQ ALAKQIHAVR YLECSALQQD GVKEVFAEAV RAVLNPTPIK RGRSC

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human RhoG protein under reducing conditions stained with Coomassie blue