

Human NCE2/UBE2F protein (Recombinant) (His-Tag) (STJP017716)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 30% Glycerol, 1mM DTT, 200mM NaCl
Molecular Weight	23.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	140739
Gene Symbol	UBE2F
UniProt ID	UBE2F_HUMAN
Immunogen Region	1-185aa
Immunogen Sequence	MLTLAS KLKRDDGLKG SRTAATASDS TRRVSVRDKL LVKEVAEELEA NLPCTCKVHF PDPNKLHCFQ LTVTPDEGY Y QGGKFQFETE VPDAYNMVPP KVKCLTKIWH PNITETGEIC LSLREHSID GTGWAPTRTL KDVVWGLNSL FTDLLNFDDP LNIEAAEHL RDKEDFRNKV DDYIKRYAR

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human NCE2/UBE2F protein under reducing conditions