

Human Ubc6/UBE2J2 protein (Recombinant) (His-Tag) (STJP017717)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 30% Glycerol, 1mM DTT, 200mM NaCl
Molecular Weight	28kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	118424
Gene Symbol	UBE2J2
UniProt ID	UB2J2_HUMAN
Immunogen Region	1-226aa
Immunogen Sequence	MSSTSS KRPTTATQR LKQDYLRKK DPVPYICAEP LPSNILEWHY VVRGPEMTPY EGGYYHGKLI FPREFPFKPP SIYMITPNGR FKCNTLRLCLS ITDFHPDTWN PAWSVSTILT GLLSFMVEKG PTLGSIETSD FTKRQLAVQS LAFNLKDKVF CELFPEVVVE IKQKQKAQDE LSSRPQTLPL PDVVPDGETH LVQNGI

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human Ubc6 UBE2J2 protein 3µg under reducing conditions with Coomassie blue stain