

Human DIMIT1 protein (Recombinant) (His-Tag) (STJP017731)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 30% Glycerol, 2mM DTT, 200mM NaCl, 2mM EDTA
Molecular Weight	37.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27292
Gene Symbol	DIMIT1
UniProt ID	DIM1_HUMAN
Immunogen Region	1-313aa
Immunogen Sequence	MPKVKSGAI GRRRGRQEQR RELKSAGGLM FNTGIGQHIL KNPLIINSII DKAALRPTDV VLEVGP GTGN MTVKLLEKAK KVVACELDPR LVAELHKRVQ GTPVASKLQV LVGDVLKTDL PFFDTCVANL PYQISSPFVF KLLHRPFFR CAILMFQREF ALRLVAKPGD KLYCRLSINT QLLARVDHLM KVGKNNFRPP PKVESS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel showing recombinant human DIMIT1 protein 3µg under reducing conditions with Coomassie blue stain