

Human Crystallin gamma N/CRYGN protein (Recombinant) (His-Tag) (STJP017736)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% Glycerol
Molecular Weight	23.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	155051
Gene Symbol	CRYGN
UniProt ID	CRGN_HUMAN
Immunogen Region	1-182aa
Immunogen Sequence	MGSSHHHHH SSSLVPRGSH MGSMAQRSG KITLYEGKHF TGQKLEVFGD CDFQDRGFM NRVNSIHVES GAWVCFNHPD FRGQQFILEH GDYPDFFRW N SHSDHMGSCR PVGMHGEHFR LEIFEGCNFT GQCLEFLEDS PFLQSRGWVK NCVNTIKVYG DGAAWSPRSF GAEDFQLSSS LQSDQGPEEA TTKPATTQPP FLTANL

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human Crystallin gamma N CRYGN protein under reducing conditions