

## Human ADPRH protein (Recombinant) (His-Tag) (STJP017742)

### GENERAL INFORMATION

|               |          |
|---------------|----------|
| Product Type  | Proteins |
| Applications  | SDS-PAGE |
| Host / Source | E.coli   |

### PRODUCT PROPERTIES


|                     |                                                                                                                                                     |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| Concentration       | 0.5 mg/mL                                                                                                                                           |
| Formulation         | Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 1mM DTT, 100mM NaCl                                                                |
| Molecular Weight    | 42.1 kDa                                                                                                                                            |
| Storage Instruction | For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                    |                                                                                                                                                                                                                                                                  |
|--------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">141</a>                                                                                                                                                                                                                                              |
| Gene Symbol        | <a href="#">ADPRH</a>                                                                                                                                                                                                                                            |
| UniProt ID         | <a href="#">ADPRH_HUMAN</a>                                                                                                                                                                                                                                      |
| Immunogen Region   | 1-357aa                                                                                                                                                                                                                                                          |
| Immunogen Sequence | MGSSHHHHHH SSGLVPRGSH MGSHEKYVA AMVLSAAGDA LGYYNGKWEF LQDGEKIHRQ<br>LAQLGGDLAL DVGRWRVSD TVMHLATAEA LVEAGKAPKL TQLYYLLAKH YQDCMEDMDG<br>RAPGGASVHN AMQLKPGKPN GWRIPFNSHE GGCGAAMRAM CIGLRFPHHS QLDTLIQVSI ESGRMTHHHP<br>TGYLALASA LFTAYAVNSR PPLQWGKGLM ELLPEAKK |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|



SDS-PAGE gel of recombinant human ADPRH protein under reducing conditions Coomassie blue stain