

Human DNAJC12 protein (Recombinant) (His-Tag) (STJP017745)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 2mM DTT, 100mM NaCl
Molecular Weight	26kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	56521
Gene Symbol	DNAJC12
UniProt ID	DJC12_HUMAN
Immunogen Region	1-198aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSMDAILN YRSEDTEYY TLLGCDELSS VEQILAEFKV RALECHPDKH PENPKAVETF QKLQKAKEIL TNEESRARYD HWRRSQMSMP FQQWEALNDS VKTSMHWVVR GKKDLMLEES DKTHTTKMEI EECNEQREK KEELASTAEK TEQKEPKPLE KSVSPQNSDS SGFADVNGWH LRFRWSDAP SELLRKFRNY EI

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human DNAJC12 protein under reducing conditions stained with Coomassie blue