

Human DBNDD1 protein (Recombinant) (His-Tag) (STJP017752)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 0.1M NaCl, 1mM DTT
Molecular Weight	19.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79007
Gene Symbol	DBNDD1
UniProt ID	DBND1_HUMAN
Immunogen Region	1-158aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSHPMEPEG AGTGEIVKEA EVPQAALGVP AAGTGDNGHT PVVEEVGGIP VPAPGLLQVT ERRQPLSSVS SLEVHFDLLD LTELTDMSDQ ELAEVFADSD DENLNTESPA GLHPLPRAGY LRSPSWTRTR AEQSHEKQPL GDPERQATVL DFTLTVERPQ ED

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human DBNDD1 protein under reducing conditions