

## Human ASF1B protein (Recombinant) (His-Tag) (STJP017753)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

### PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol 0.15M NaCl, 1mM DTT
Molecular Weight	23.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">55723</a>
Gene Symbol	<a href="#">ASF1B</a>
UniProt ID	<a href="#">ASF1B_HUMAN</a>
Immunogen Region	1-202aa
Immunogen Sequence	MAKVSVLNVA VLENPSPFHS PFRFEISFEC SEALADDLEW KIIYVGSAES EEFQILDVSV LVGPVPAGRH MFVFAQADAPN PSLIPETDAV GVTVVLTICT YHGQEFIRVG YYVNNYLNPELRENPPMKP DFSQLQRNIL ASNPRVTRFH INWDNNMDRL EAIETQDPSL GCGLPLNCTP IKGLGLPGCI PGLLPENSMD CILE

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human ASF1B protein 3µg under reducing conditions Coomassie blue stain