

Human JNK2/MAPK9 protein (Recombinant) (His-Tag) (STJP017758)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 10% Glycerol, 1mM DTT
Molecular Weight	46.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5601
Gene Symbol	MAPK9
UniProt ID	MK09_HUMAN
Immunogen Region	1-382aa
Immunogen Sequence	MSDSKC DSQFYSVQVA DSTFTVLKRY QQLKPIGSGA QGIVCAAFDT VLGINVAVKK LSRPFQNTTH AKRAYRELVL LKCVNHKNII SLLNVFTPQK TLEEFQDVYL VMELMDANLC QVIHMELDHE RMSYLLYQML CGIKHLHSAG IHRDLKPSN IVKSDCTLK ILDFGLARTA CTNFMMPYV VTRYRAPEV ILGMGY

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing 3µg recombinant human JNK2 MAPK9 protein under reducing conditions stained with Coomassie blue