

Human Ephrin-B3 protein (Recombinant) (His-Tag) (STJP017767)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol, 0.1M NaCl, 2M urea
Molecular Weight	24.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1949
Gene Symbol	EFNB3
UniProt ID	EFNB3_HUMAN
Immunogen Region	28-226aa
Immunogen Sequence	LSLEP VYWSANKRF QAEGGYVLYP QIGDRLDLLC PRARPPGPHS SPNYEFYKLY LVGGAQGRRC EAPPAPNLLL TCDRPDLRLR FTIKFQEYSP NLWGHEFRSH HDYYIIATSD GTREGLESLO GGVCLTRGMK VLLRVGQSPR GGAVPRKPVs EMPMERDRGA AHSLEPGKEN LPGDPTSNAT SRGAEGPLPP PSMP

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human Ephrin-B3 protein under reducing conditions stained with Coomassie blue