

Human 2B4/CD244 protein (Recombinant) (His-Tag) (STJP017772)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea
Molecular Weight	25.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51744
Gene Symbol	CD244
UniProt ID	CD244_HUMAN
Immunogen Region	19-224aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSHGKGCQG SADHVVSISG VPLQLQPNSI QTKVDSIAWK KLLPSQNGFH HILKWENGSL PSNTSND RFS FIVKNLSLLI KAAQQQDSGL YCLEVTSISG KVQTATFQVF VFDPKVEKPRL QGQGGKILDRG RCQVALSCLV SRDGNVSYAW YRGSKLIQTA GNLTYLDEEV DINGTHYTC NVSNPVSWEH HTLNLTQDCQ NAHQEFRF

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human 2B4/CD244 protein under reducing conditions