

Human C4BPB protein (Recombinant) (His-Tag) (STJP017773)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 0.15M NaCl
Molecular Weight	29kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	725
Gene Symbol	C4BPB
UniProt ID	C4BPB_HUMAN
Immunogen Region	18-252aa
Immunogen Sequence	SDAEH CPELPPVDNS IFVAKEVEGQ ILGTYVCIKG YHLVGKKTLF CNASKEWDNT TTECR LGHCP DPVLVNGEFS SSGPVNVSDK ITFMCNDHYI LKGSNRSQCL EDHTWAPPPF ICKSRDCDPP GNPVHG YFEG NNFTLGSTIS YYCEDRYLV GVQEQQCV DG EWSSALPVCK LIQEAPKPEC EKALLAFQES KNLCEA

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human C4BPB protein under reducing conditions stained with Coomassie blue