

Human ARH3/ADPRS protein (Recombinant) (His-Tag) (STJP017774)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 30% Glycerol, 1mM DTT
Molecular Weight	41.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54936
Gene Symbol	ADPRS
UniProt ID	ADPRS_HUMAN
Immunogen Region	1-363aa
Immunogen Sequence	MAAAAM AAAAGGGAGA ARSLSRFRGC LAGALLGDCV GSFYEAHDTV DLTSVLRHVQ SLEPDGTPG SERTEALYYT DDTAMARALV QSLLAKEAFD EVDMAHRFAQ EYKDPDRGY GAGVVTVFKK LLNPKCRDVF EPARAFNGK GSYGNGGAMR VAGISLAYSS VQDVQKFARL SAQLTHASSL GYNGAILQAL AVHLAL

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human ARH3/ADPRS protein under reducing conditions
Coomassie blue stain