

## Human PAFAH1B3 protein (Recombinant) (His-Tag) (STJP017783)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

### PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 10% Glycerol
Molecular Weight	28.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">5050</a>
Gene Symbol	<a href="#">PAFAH1B3</a>
UniProt ID	<a href="#">PA1B3_HUMAN</a>
Immunogen Region	1-231aa
Immunogen Sequence	MSGEENP ASKPTPVQDV QGDGRWMSLH HRFVADSKDK EPEVFIGDS LVQLMHQCEI WRELFSP LHA LNF GIGGDGT QHVLWRENG ELEHIRPKIV VVWVG TNNHG HTAEQVTGGI KAI VQLVNER QPQARVVVLG LLPRGQHPNP LREKNRQVNE LVRAALAGHP RAHFLDADPG FVHSDGTISH HDMYDYLHLS RLG YTP

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



STJP017783 recombinant human PAFAH1B3 protein SDS-PAGE under reducing conditions with Coomassie blue staining