

## Human Coagulation Factor III/Tissue Factor protein (Recombinant) (His-Tag) (STJP017792)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

### PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 10% Glycerol, 1mM DTT
Molecular Weight	25.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">2152</a>
Gene Symbol	<a href="#">F3</a>
UniProt ID	<a href="#">TF_HUMAN</a>
Immunogen Region	33-251aa
Immunogen Sequence	MSGTTNTVAA YNLTKWSTNF KTILEWEPKP VNQVYTVQIS TKSGDWKSKC FYTTDTECDL TDEIVKDVKQ TYLARVFSYP AGNVESTGSA GEPLYENSPE FTPYLETNLG QPTIQSFEQV GTKVNVTVED ERTLVRNNNT FLRLRDVFGK DLIYTLYYWK SSSSGKKTAK TNTNEFLIDV DKGENYCFVS QAVIPSRTVN RKSTDSPVEC MGQEKGEFRE

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human Coagulation Factor III Tissue Factor protein under reducing conditions