

Human DKK2 protein (Recombinant) (His-Tag) (STJP017809)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% Glycerol
Molecular Weight	27.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27123
Gene Symbol	DKK2
UniProt ID	DKK2_HUMAN
Immunogen Region	34-259aa
Immunogen Sequence	KLNSIK SSLGGETPGQ AANRSAGMYQ GLAFGGSKKG KNLGQAYPCS SDKECEVGRY CHSPHQGSSA CMVCRKKKKR CHRDMCCPS TRCNNGICIP VTESILTPHI PALDGTRHRD RNHGHYSNHD LGWQNLGRPH TKMSHIKGHE GDPCLRSSDC IEGFCCARHF WTKICKPVLH QGEVCTKQRK KGSHGLEIFQ RCDDAK

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE of recombinant human DKK2 protein under reducing conditions Coomassie blue stain