

## Human SEC22B protein (Recombinant) (His-Tag) (STJP017818)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

### PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 2M urea, 10% Glycerol, 0.15M NaCl
Molecular Weight	23.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">9554</a>
Gene Symbol	<a href="#">SEC22B</a>
UniProt ID	<a href="#">SC22B_HUMAN</a>
Immunogen Region	14-194aa
Immunogen Sequence	LPLAAS MQEDEQSGRD LQQYQSQAQK LFRKLNEQSP TRCTLEAGAM TFHYIIEQGV CYLVLCEAAF PKKLAFAYLE DLHSEFDEQH GKKVPTVSRP YSFIEFDTFI QKTKKLYIDS RARRNLGSIN TELQDVQRIM VANIEEVLQR GEALSALDSK ANNLSLSSKK YRQDAKYLMN RSTYA

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human SEC22B protein under reducing conditions with Coomassie blue stain