

Human MOCS2 protein (Recombinant) (His-Tag) (STJP017821)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 50% Glycerol, 1mM DTT
Molecular Weight	25 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4338
Gene Symbol	MOCS2
UniProt ID	MOC2B_HUMAN
Immunogen Region	1-188aa
Immunogen Sequence	MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSMSL EISSSCFSLE TKLPLSPPLV EDSAFEPSPRK DMDEVEEKSK DVINFTAEKL SVDEVSQLVI SPLCGAISLF VGTTRNNFEG KKVISLEYEA YLPMAENEVR KICSDIRQKW PVKHIAVFHR LGLVPVSEAS IIIAVSSAHR AASLEAVSYA IDTLKAKVPI WKKEIYEES TWKGNKECFW ASNS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human MOCS2 protein under reducing conditions stained with Coomassie blue