

Human SH3BGRL3 protein (Recombinant) (His-Tag) (STJP017824)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% Glycerol, 1mM DTT
Molecular Weight	12.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	83442
Gene Symbol	SH3BGRL3
UniProt ID	SH3L3_HUMAN
Immunogen Region	1-93aa
Immunogen Sequence	MSGLRVY STSVTGSREI KSQQSEVTRI LDGKRIQYQL VDISQDNALR DEMRALAGNP KATPPQIVNG DQYCGDYELF VEAVEQNTLQ EFLKLA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human SH3BGRL3 protein Coomassie blue stain reducing conditions