

Human ANKRD1 protein (Recombinant) (His-Tag) (STJP017827)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.5) containing 0.2M NaCl, 50% Glycerol, 5mM DTT
Molecular Weight	38.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27063
Gene Symbol	ANKRD1
UniProt ID	ANKR1_HUMAN
Immunogen Region	1-319aa
Immunogen Sequence	MVLKVE ELVTGKKNGN GEAGEFLPED FRDGEYEAAV TLEKQEDLKT LLAHPVTLGE QQWKSEKQRE AELKKKKLEQ RSKLENLEDL EIIQLKKRK KYRKTQVPV KEPEPEIITE PVDVPTFLKA ALENKLPVVE KFLSDKNNPD VCDEYKRTAL HRACLEGH LA IVEKLMEAGA QIEFRDMLES TAIHWASRGG NLDVLK

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human ANKRD1 protein 3µg under reducing conditions Coomassie blue stain