

Human p21/CIP1/CDKN1A protein (Recombinant) (His-Tag) (STJP017834)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 2M urea, 10% Glycerol, 0.1M NaCl
Molecular Weight	20.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1026
Gene Symbol	CDKN1A
UniProt ID	CDN1A_HUMAN
Immunogen Region	1-164aa
Immunogen Sequence	MSEPAGDVRQ NPCGSKACRR LFGPVDSEQL SRDCDALMAG CIQEAREERWN FDFVTETPLE GDFAWERVRG LGLPKLYLPT GPRRGRDELG GGRRPGTSPA LLQGTAEEDH VDLSSLCTLV PRSGEQAEGS PGGPGDSQGR KRRQTSMTDF YHSKRRLIFS KRKP

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human p21 CIP1 CDKN1A protein under reducing conditions