

Human GNAZ protein (Recombinant) (His-Tag) (STJP017836)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% Glycerol
Molecular Weight	43.0 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2781
Gene Symbol	GNAZ
UniProt ID	GNAZ_HUMAN
Immunogen Region	1-355aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGCQRSSEK EAARRSRID RHLRSESQRQ RREIKLLLG TSNSGKSTIV KQMKIIHSGG FNLEACEYK PLIYNAIDS LTRIIRALAA LRIDFHNPDY AYDAVQLFAL TGPAESKGEI TPELLGVMRR LWADPGAQAC FSRSEYHLE DNAAYYLDL ERIAADYIP TVEDILRSRD MTTGIVENKF TFKELTFKMV DVGQQRSE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human GNAZ protein 3µg under reducing conditions
Coomassie blue stain