

Human COMMD9 protein (Recombinant) (His-Tag) (STJP017856)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% Glycerol, 1mM DTT
Molecular Weight	24.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29099
Gene Symbol	COMMD9
UniProt ID	COMD9_HUMAN
Immunogen Region	1-198aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSMAALTAE HFAALQSLK ASSKDVVRQL CQESFSSAL GLKLLDVTC SSLSVTQEEA EELLQALHRL TRLVAFRDLS SAEAILALFP ENFHQNLKLN LTKIILEHVS TWRTEAQANQ ISLPRLVDLD WRVDIKTSSD SISRMVPTC LLQMKIQEDP SLCGDKPSIS AVTVELSKET LDTMLDGLGR IRDQLSAVAS K

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human COMMD9 protein under reducing conditions stained with Coomassie blue