

Human NDUFA2 protein (Recombinant) (His-Tag) (STJP017878)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% Glycerol, 1mM DTT
Molecular Weight	13.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4695
Gene Symbol	NDUFA2
UniProt ID	NDUA2_HUMAN
Immunogen Region	1-99aa
Immunogen Sequence	MAAAAAS RGVGAKLGLR EIRIHLQQRS PGSQGVRDFI EKRYVELKKA NPDLPIIRE CSDVQPKLWA RYAFGQETNV PLNNFSADQV TRALENVLSG KA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE of recombinant human NDUFA2 protein under reducing conditions stained with Coomassie blue