

Human RNA polymerase II subunit RPB11-a/POLR2J protein (Recombinant) (His-Tag) (STJP017882)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 40% Glycerol, 5mM DTT, and 2mM EDTA
Molecular Weight	15.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5439
Gene Symbol	POLR2J
UniProt ID	RPB11_HUMAN
Immunogen Region	1-117aa
Immunogen Sequence	MNAPPAF ESFLFEGEK KITINKDTKV PNACFLTINK EDHTLGNIHK SLLKDPQVL FAGYKVPHPLEHKIIIRVQT TPDYSPQEAFTNAITDLISE LSLLEERFRV AIKDKQEGIE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human RNA polymerase II subunit RPB11-a POLR2J protein