

Human NDUFAB1 protein (Recombinant) (His-Tag) (STJP017883)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 40% Glycerol, 1mM DTT, 2mM EDTA
Molecular Weight	37 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51103
Gene Symbol	NDUFAB1
UniProt ID	CIA30_HUMAN
Immunogen Region	25-327aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSYPFLGIR FAEYSSSLQK PVASPGKASS QRKTEGDLQG DHQKEVALDI TSSEKPDVS FDKAIRDEAI YHFRLLKDEI VDHWRGPEGH PLHEVLLEQA KVVWQFRGKE DLDKWTVTSD KTIGGRSEVF LKMGKNNQSA LLYGTLSSA PQDGESTRSG YCAMISRIPR GA FERKMSYD WSQFNTLYLR VRGDGRPWMV NIKEDTDF

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human NDUFAB1 protein under reducing conditions with Coomassie blue stain