

Human PDZD11 protein (Recombinant) (His-Tag) (STJP017887)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% Glycerol, 1mM DTT
Molecular Weight	18.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51248
Gene Symbol	PDZD11
UniProt ID	PDZ11_HUMAN
Immunogen Region	1-140aa
Immunogen Sequence	MGSMDSRIPY DDYPVFLPA YENPPAWIPP HERVHHPDYN NELTQFLPRT ITLKPPGAQ LGFNIRGGKA SQLGIFISKV IPDSDAHRAG LQEGDQVLAV NDVDFQDIEH SKAVEILKTA REISMRVRFY PNYHRQKER TVH

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human PDZD11 protein under reducing conditions Coomassie blue stain