


Human SLD5/GINS4 protein (Recombinant) (His-Tag) (STJP017900)

GENERAL INFORMATION	
Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli
PRODUCT PROPERTIES	
Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 20% Glycerol, 1mM DTT
Molecular Weight	28.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	84296
Gene Symbol	GINS4
UniProt ID	SLD5_HUMAN
Immunogen Region	1-223aa
Immunogen Sequence	MGSSHHHHH SSGLVPRGSH MTEEVDLFGQ DSDGGSEEVV LTPAELIERL EQAWMNEKFA PELLESKPEI VECVMEQLEH MEENLRRAKR EDLKVSIHQM EMERIRYVLS SYLRCLMKI EKFFPHVLEK EKTRPEGEPS SLSPLELAFA REFMANATESY LKNVALKHMP PNLQKVDLFR AVPKPDLSY VFLRVRRERQE NILVEPDTDE QRDYVIDLEK GSQHLIRY
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human SLD5/GINS4 protein under reducing conditions
Coomassie blue stain