

Human Endophilin B2/SH3GLB2 protein (Recombinant) (His-Tag) (STJP017922)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 30% Glycerol, 1mM DTT
Molecular Weight	46.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	56904
Gene Symbol	SH3GLB2
UniProt ID	SHLB2_HUMAN
Immunogen Region	1-395aa
Immunogen Sequence	MDFNMKK LASDAGIFFT RAVQFTEEFK GQAEKTELDA HFENLLARAD STKNWTEKIL RQTEVLLQPN PSARVEEFY EKLDKVPSPR VTNGELLAQY MADAASELGP TTPYGKTLIK VAAAEKQLGA AERDFIHTAS ISFLTPLRNF LEGDWKTISK ERRLLQNRRL DLDACKARLK KAKAAEAKAT TVPDFQETRP RNYILS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human Endophilin B2 SH3GLB2 protein under reducing conditions