

## Human RPS5 protein (Recombinant) (His-Tag) (STJP017934)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

### PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 50% Glycerol, 2mM DTT
Molecular Weight	25.3 Da
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">6193</a>
Gene Symbol	<a href="#">RPS5</a>
UniProt ID	<a href="#">RS5_HUMAN</a>
Immunogen Region	1-204aa
Immunogen Sequence	MTEWETA APAVAETPDI KLFQKWTDD VQINDISLQD YIAVKEYAK YLPHSAGRYA AKRFRKAQCP IVERLTNSMM MHGRNNGKKL MTRIVKHAF EIIHLLTGEN PLQVLVNII NSGPREDSTR IGRAGTVRRQ AVDVSPLRRV NQAIWLLCTG AREAAFRNIK TIAECLADEL INAAKGSSNS YAIKKKDELE RVAKSN

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human RPS5 protein 3µg under reducing conditions stained with Coomassie blue