

Human RNA polymerase II subunit RPB9/POLR2I protein (Recombinant) (His-Tag) (STJP017935)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 20% Glycerol, 1mM DTT
Molecular Weight	17 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5438
Gene Symbol	POLR2I
UniProt ID	RPB9_HUMAN
Immunogen Region	1-125aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSMEPDGTY EPGFVGIRFC QECNNMLYPK EDKENRILLY ACRNCDYQQE ADNSCIYVVK ITHEVDELTDQ IADVSQDPT LPRTEHDPCQ KCGHKEAVFF QSHSARAEDA MRLYYVCTAP HCGHRWTE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human RNA polymerase II subunit RPB9 POLR2I protein under reducing conditions