

Human IRGM protein (Recombinant) (His-Tag) (STJP017940)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 2M urea, 10% Glycerol
Molecular Weight	20.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	345611
Gene Symbol	IRGM
UniProt ID	IRGM_HUMAN
Immunogen Region	23-181aa
Immunogen Sequence	KETLKIVSR TPNITMAGD SGNMGSTFIS ALRNTGHEGK ASPPTLVKA TQRCASYFSS HFSNVVLWDL PGTGSATTL ENYLMEMQFN RYDFIMVASA QFSMNHVMLA KTAEDMGKKF YIYWTKLDMD LSTGALPEVQ LLQIRENVLE NLQKERVCEY

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human IRGM protein under reducing conditions Coomassie blue stain