

Human ICT1 protein (Recombinant) (His-Tag) (STJP017941)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 30% Glycerol, 1mM DTT
Molecular Weight	22.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3396
Gene Symbol	MRPL58
UniProt ID	ICT1_HUMAN
Immunogen Region	30-206aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGS LHKQKDG TEFKSIYSLD KLYPESQGSD TAWRVPNGAK QADSDIPLDR LTISYCRSSG PGGQNVNKNV SKAEVRFHLA TAEWIAEPVR QKIAITHKNK INRLGELILT SESSRYQFRN LADCLQKIRD MITEASQTPK EPTKEDVKLH RIRIENMNRE RLRQKRIHSA VKTSRRVDM

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human ICT1 protein under reducing conditions stained with Coomassie blue