

## Human GLTP protein (Recombinant) (His-Tag) (STJP017956)

### GENERAL INFORMATION

|               |          |
|---------------|----------|
| Product Type  | Proteins |
| Applications  | SDS-PAGE |
| Host / Source | E.coli   |

### PRODUCT PROPERTIES


|                     |                                                                                                                                                     |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| Concentration       | 0.5 mg/mL                                                                                                                                           |
| Formulation         | Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 20% Glycerol, 1mM DTT                                                                 |
| Molecular Weight    | 26.2 kDa                                                                                                                                            |
| Storage Instruction | For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                    |                                                                                                                                                                                                                                                                    |
|--------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">51228</a>                                                                                                                                                                                                                                              |
| Gene Symbol        | <a href="#">GLTP</a>                                                                                                                                                                                                                                               |
| UniProt ID         | <a href="#">GLTP_HUMAN</a>                                                                                                                                                                                                                                         |
| Immunogen Region   | 1-209aa                                                                                                                                                                                                                                                            |
| Immunogen Sequence | MGSSHHHHH SSGLVPRGSH MGSMALLAEH LLKPLPADKQ IETGPFLEAV SHLPPFFDCL GSPVFTPIKA<br>DISGNITKIK AVYDTNPAKF RTLQNIILEVE KEMYGAEWPK VGATLALMWL KRGLRFIQVF LQSIDGGERD<br>ENHPNLIRVN ATKAYEMALK KYHGWIVQKI FQAALYAAPY KSDFLKALSK GQNVTEEECL EKIRLFLVNY<br>TATIDVIYEM YTMNAEL |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|



SDS-PAGE gel of recombinant human GLTP protein 3µg under reducing conditions Coomassie blue stain