

Human IFN-alpha 7/IFNA7 protein (Recombinant) (His-Tag) (STJP017980)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% Glycerol
Molecular Weight	22.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3444
Gene Symbol	IFNA7
UniProt ID	IFNA7_HUMAN
Immunogen Region	24-189aa
Immunogen Sequence	MCDLPQ THSLRNRRL ILLAQMGRI PFCLKDRHE FRFPEEEFDG HQFQKTQAIS VLHEMIQQTFLFSTEDSSA AWEQSLEKF STELYQQLND LEACVIQEVG VEETPLMNED FILAVRKYFQ RITLYLMEKK YSPCAWEVVR AEIMRSFSFS TNLKKGLRRK D

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human IFN-alpha 7 IFNA7 protein under reducing conditions