

Human JAM-B/VE-JAM protein (Recombinant) (His-Tag) (STJP017981)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% Glycerol
Molecular Weight	26.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	58494
Gene Symbol	JAM2
UniProt ID	JAM2_HUMAN
Immunogen Region	21-238aa
Immunogen Sequence	MGSSHHHHH SSGLVPRGSH MGSLGYHKAY GFSAPKDQQV VTAVEYQEI LACKTPKKT V SSRLEWKKLG RSVSFVYYQQ TLQGDFKNRA EMIDFNIRIK NVTRSDAGKY RCEVSAPSEQ GQNLEEDTVT LEVLVAPAVP SCEVPSSALS GTVVELRCQD KEGNPAPEYT WFKDGIRLLE NPRLGSQSTN SSYTMNTKTG TLQFNTVSKL DTGEYSCEAR NSVGYYRC

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel showing 3µg recombinant human JAM-B/VE-JAM protein under reducing conditions