

Human XRCC3 protein (Recombinant) (His-Tag) (STJP017983)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol 0.4M urea
Molecular Weight	40 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7517
Gene Symbol	XRCC3
UniProt ID	XRCC3_HUMAN
Immunogen Region	1-346aa
Immunogen Sequence	MDLDLLDLNP RIIAAIKKAK LKSVKEVLHF SGPDLKRLTN LSSPEVWHLL RTASLHLRGS SILTALQLHQ QKERFPTQHQ RLSLGCPVLD ALLRGGLPLD GITELAGRSS AGKTQLALQL CLAVQFFPRQH GGLEAGAVYI CTEDAFPHKR LQQLMAQQPR LRTDVP GELL QKLRFQSQIF IEHVADVDTL LECVNKKVPV LLSRGM

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human XRCC3 protein under reducing conditions Coomassie blue stain