

Human Fructosamine-3-kinase/FN3K protein (Recombinant) (His-Tag) (STJP017991)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% Glycerol, 1mM DTT
Molecular Weight	37 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	64122
Gene Symbol	FN3K
UniProt ID	FN3K_HUMAN
Immunogen Region	1-309aa
Immunogen Sequence	MEQLLRA ELRTATLRAF GPGGAGCISE GRAYDTDAGP VFKVNRRTQ ARQMFEGEVA SLEALRSTGL VRVPRPMKVI DLPGGGAFAV MEHLKMKSL SQASKLGEQM ADLHLYNQKL REKLKEEENT VGRRGEGAEP QYVDKFGFHT VTCCGFIPQV NEWQDDWPTF FARHRLQAQL DLIEKDYADR EARELWSRLQ VKIPDL

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human Fructosamine-3-kinase FN3K protein under reducing conditions