

Human FN3KRP protein (Recombinant) (His-Tag) (STJP017999)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% Glycerol
Molecular Weight	36.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79672
Gene Symbol	FN3KRP
UniProt ID	KT3K_HUMAN
Immunogen Region	1-309aa
Immunogen Sequence	MEELLRR ELGCSSVRAT GHSGGGCISQ GRSYDTDQGR VFKVNPKE ARRMEFEGEMA SLTAILKTNT VKVPKPIKVL DAPGGGSLV MEHMDMRHLS SHAAKLGAQL ADLHLDNKKL GEMRLKEAGT VGRGGGQEE PFVARFGFDV VTCCGYLPQV NDWQEDWVVF YARQRIQPQM DMVEKESGDR EALQLWSALQ LKIPDL

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE of recombinant human FN3KRP protein 3µg under reducing conditions Coomassie blue stain