

Human Inhibin A protein (Recombinant) (His-Tag) (STJP018005)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% Glycerol
Molecular Weight	17 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3623
Gene Symbol	INHA
UniProt ID	INHA_HUMAN
Immunogen Region	233-366aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MSTPLMSWPW SPSALRLLQR PPEEPAAHAN CHRVALNISF QELGWERWIV YPPSFIFHYC HGGCGLHIPP NLSLPVPGAP PTPAQPYSL PGAQPCCAAL PGTMRPLHVR TTSDGGYSFK YETVPNLLTQ HCACI

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human Inhibin A protein 3 µg under reducing conditions stained with Coomassie blue