

Human GALT protein (Recombinant) (His-Tag) (STJP018011)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.5) containing 0.2M NaCl, 10% Glycerol
Molecular Weight	45.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2592
Gene Symbol	GALT
UniProt ID	GALT_HUMAN
Immunogen Region	1-379aa
Immunogen Sequence	MSRSGT DPQQRQQASE ADAAAATFRA NDHQHIRYNP LQDEWVLVSA HRMKRPWQQGQ VEPQLLKTVP RHDPLNPLCP GAIRANGEVN PQYDSTFLFD NDFPALQPD A PSPGPSDHPL FOAKSARGVC KVMCFHPWSD VTLPLMSVPE IRAVVD AWAS VTEELGAQYP WVQIFENKGA MMGCSNPHPH CQVWASSFLP DIAQRE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human GALT protein 3µg under reducing conditions stained with Coomassie blue