

## Human DNAJB4 protein (Recombinant) (His-Tag) (STJP018019)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

### PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 7.5) containing 0.2M NaCl, 30% Glycerol, 1mM DTT
Molecular Weight	40 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">11080</a>
Gene Symbol	<a href="#">DNAJB4</a>
UniProt ID	<a href="#">DNJB4_HUMAN</a>
Immunogen Region	1-337aa
Immunogen Sequence	MGKDYYC ILGIEKGASD EDIKKAYRKQ ALKFHPDKNK SPQAEKFK E VAEAYEVLSD PKKREIYDQF GEEGLKGGAG GTDGQGGTFR YTFHGDPHAT FAAFFGGSNP FEIFFGRRMG GGRDSEEMEI DGDPFSAFGF SMNGYPRDRN SVGPSRLKQD PPVIHELKRV LEEIYSGCTK RMKISRKRLN ADGRSYRSED KILTIE

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human DNAJB4 protein 3µg under reducing conditions with Coomassie blue stain