

Human OBFC1 protein (Recombinant) (His-Tag) (STJP018042)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% Glycerol, 1mM DTT
Molecular Weight	44.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79991
Gene Symbol	STN1
UniProt ID	STN1_HUMAN
Immunogen Region	1-368aa
Immunogen Sequence	MQPGSSR CEEETPSLLW GLDPVFLAFA KLYIRDILDM KESRQVPGVF LYNGHPIKQV DVLGTVIGVR ERDAFYSGV DDSTGVINCI CWKKLNTESV SAAPSAAREL SLTSQKKLQ ETIEQKTKIE IGDITRVRS IRTYREEREI HATAYYKVD PVWNIQIARM LELPTIYRKV YDQPFHSSAL EKEEALSNPG ALDLPS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human OBFC1 protein under reducing conditions Coomassie blue stain