

Human GOT2 protein (Recombinant) (His-Tag) (STJP018050)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% Glycerol, 1mM DTT
Molecular Weight	47 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2806
Gene Symbol	GOT2
UniProt ID	AATM_HUMAN
Immunogen Region	30-430aa
Immunogen Sequence	SSWWTHV EMGPPDPILG VTEAFKRDTN SKKMNLGVGA YRDDNGKPYV LPSVRKAEAQ IAAKNLDKEY LPIGGLAEFC KASAEALALGE NSEVLKSGRF VTVQTISGTG ALRIGASFLQ RFFKFSDRVF LPKPTWGNHT PIFRDAGMQL QGYRYYPKT CGFDFTGAVE DISKIQSV LLLHACAHNP TGVDPRPEQW KEIATV

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human GOT2 protein under reducing conditions stained with Coomassie blue