

## Human RNA polymerase II subunit RPB11-b1/POLR2J2 protein (Recombinant) (His-Tag) (STJP018059)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

### PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% Glycerol, 1mM DTT, 250mM Imidazole
Molecular Weight	15.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">246721</a>
Gene Symbol	<a href="#">POLR2J2</a>
UniProt ID	<a href="#">RPB1B_HUMAN</a>
Immunogen Region	1-115aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGS MNAPP AF ESFL FEGEK ITINKDTKVP KACLFTINKE DHTLGNI IKS QLLKDPQVLF AGYKVPHPLE HKIIIRVQTT PDYSPQE AFT NAITDLISEL SLLEERFRTC LLPLRLLP

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human RNA polymerase II subunit RPB11-b1 POLR2J2 protein