

Human AMSH-LP protein (Recombinant) (His-Tag) (STJP018062)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 20% Glycerol, 0.1mM PMSF
Molecular Weight	52.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57559
Gene Symbol	STAMBPL1
UniProt ID	STALP_HUMAN
Immunogen Region	1-436aa
Immunogen Sequence	MDQPFTVN SLKKLAAMPD HTDVLSPEE RVRALSKLGC NITISEDITP RRYFRSGVEM ERMASVYLEE GNLENAFVLY NKFITLFVEK LPNHRDYQQC AVPEKQDIMK KLKEIAFPRT DELKNDLLKK YNVEYQEYLQ SKNKYKAEIL KKLEHQRLIE AERKRIAQMR QQLESEQFL FFEDQLKKQE LARGQMRSSQ TSGLSE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human AMSH-LP protein under reducing conditions stained with Coomassie blue