

Human Cytokeratin 20/KRT20 protein (Recombinant) (His-Tag) (STJP018067)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% Glycerol
Molecular Weight	50.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54474
Gene Symbol	KRT20
UniProt ID	K1C20_HUMAN
Immunogen Region	1-424aa
Immunogen Sequence	MDFSRRS FHRSLSSSLQ APVVSTVGMQ RLGTTPSVYG GAGGRGIRIS NSRHTVNYGS DLTGGGDLFV GNEKMAMQNL NDRLASYLEK VRTLEQSNK LEVQIKQWYE TNAPRAGR DY SAYYRQIEEL RSQIKDAQ LQ NARCVLQIDN AKLAAEDFRL KYETERGIRL TVEADLQGLN KVFDLTLHK TDLEIQIEEL NKDLAL

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE of recombinant human Cytokeratin 20 KRT20 protein under reducing conditions